Human Gastric Tissue

Helico*Stat® Rapid Helicobacter Stain Kit Procedure

100ml Kit Item #: KTHST

Liter Kit Item#: KTHSTLT

Pint Kit Item #: KTHSTPT

Gallon Kit Item#: N/A

Control Slide(s)	Item#	Included Components	
Helicobacter Pylori	CSH0125P	Toluidine Blue Stock	Sodium Metabisulfite Solution
		Sodium Hydroxide Solution	Alcian Yellow Solution
		Periodic Acid Solution	

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify *Helicobacter pylori* bacteria. *Helicobacter pylori* are stained blue, background shades of blue, and mucin yellow.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 3-4µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

TOLUIDINE BLUE WORKING SOLUTION PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	50ml	Distilled Water	Into a chemically clean container or new/unused plasticware.
2	Add	10 drops	Toluidine Blue Stock	Mix thoroughly.
3	Add	3 drops	Sodium Hydroxide Solution	Mix thoroughly.

STAINING PROCEDURE:

Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	Immerse	Periodic Acid Solution	--	10	--	Once complete, rinse in running tap water (1 minute) and continue. Do not reuse this solution.
5	Immerse	Sodium Metabisulfite Solution	--	5	--	Once complete, rinse in running tap water (1 minute) and continue.
6	Immerse	Alcian Yellow Solution	--	5-10	--	Once complete, rinse in running tap water (1 minute) and continue.
7	Immerse	Toluidine Blue Working Solution	--	3-5	--	Once complete, rinse in running tap water (1 minute) and continue. Do not reuse this solution.
8	Air Dry	On Benchtop, or in Empty Bath	--	5-10	--	Until section is completely dry.
9	Clear	Xylene or Substitute	--	--	1	5-10 dips for 1 second each.
10	Coverslip	Permanent Mounting Media	--	--	--	

1. The Journal of Histotechnology; Vol 19, No 2; 1996, 131 - 132.
2. With modifications by AMTS R&D Department, 1979-2018.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. Standard Toluidine Blue Working Solution preparation procedure yields approximately 50ml of solution and must be scaled up to accommodate desired autostainer bath size. A minimum of 5 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 12 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	330	14	250ml Glass Stain Dish	2	248	124
30ml Glass Coplin Jar	16	317	20	200ml Bath Autostainer	2	211	106
40ml Hellendahl Jar	12	317	27	400ml Bath Autostainer	1	224	225

CE MARKINGS AND DESIGNATIONS:

	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
	Batch Code		Use By		Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use			
	Corrosive		Irritant		Health Hazard	

1. The Journal of Histotechnology; Vol 19, No 2; 1996, 131 - 132.
2. With modifications by AMTS R&D Department, 1979-2018.

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COMPOSANTS INCLUS: Toluidine Blue Stock, Sodium Hydroxide Solution, Periodic Acid Solution, Sodium Metabisulfite Solution, Alcian Yellow Solution.

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