

Accu-Edge®

Disposable Microtome Blades

Considered the best in the world
by many Histotechnicians.

Sharp, Consistent, and considered the best by many Histotechnicians, Accu-Edge® disposable microtome blades deliver excellent results in the world of microtomy. The Accu-Edge® disposable microtome blades can section down to 1 micron without striation, compression, distortion, or chattering.



Low- Profile

Pack (50 blades):
Item#: BLAEDLPPK

Case (10 packs):
Item#: BLAEDLPCS



High- Profile

Pack (50 blades):
Item#: BLAEDHPPK

Case (10 packs):
Item#: BLAEDHPCS

The dust-free blade dispenser for the Accu-Edge® Disposable microtome blades delivers microtome blades consistently in seconds and automatically positions new blades safely and precisely. These holders fit any rotary microtome, including Leica, Shandon, Leitz, Surgipath and Tissue-Tek. Each set includes a blade holder, one pack of Accu-Edge® Disposable Blades, maintenance tools, replacement parts, and a storage case.



Low-Profile Microtome Blade Holder

Maintenance Free:
Item#: BL958

Original (One-Touch):
Item#: BL962

The Blade holder for AO/Reichert Cryostat has a right-mounted plate to accommodate the off-center blade position. The Universal blade holder has a centered blade fixing plate and fits I.E.C., Bright, Leitz, Lipshaw, and other cryostat models. Both holders feature a screw-down type blade clamping, a black coated fixing plate, and an anti-roll device. Each Set includes the blade holder, one pack of the Accu-Edge low-profile disposable blades, maintenance tools, replacement parts, and a storage case.



American MasterTech
scientific laboratory supplies



Cryostat Blade Holder

Low-Profile Cryostat Blade Holder:
Item#: BL952

Universal Cryostat Blade Holder:
Item#: BL954

Copyright ©2011 American MasterTech Scientific Incorporated. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.