

DURAEEDGE®

Disposable Microtome Blades



The right decision when
you want precision.

DURAEEDGE® MICROTOME BLADES are the right decision when you want precision. The DURAEEDGE® microtome blade is manufactured to the highest standards for sharpness, consistency and durability. Our special proprietary process for hardened stainless steel ensures the quality of each finely honed and polished blade to give you a flaw-free edge. Uniform PTFE coating reduces friction, eliminates striations and compression, and requires no break-in period. Capable of sectioning tissue for biopsies in ultra-thin sections, DURAEEDGE® blades give you unsurpassed precision in the lab.



American MasterTech
scientific laboratory supplies



Low- Profile (PTFE Coated)

Pack (50 blades):
Item#: BLM00101P

Case (10 packs):
Item#: BLM00101C



High- Profile (PTFE Coated)

Pack (50 blades):
Item#: BLM00102P

Case (10 packs):
Item#: BLM00102C



High- Profile (Ceramic Coated)

Pack (50 blades):
Item#: BLM00103P

Case (10 packs):
Item#: BLM00103C

Copyright ©2011 American MasterTech Scientific Incorporated. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.