


HARRIS UNI-CORE SAMPLING TOOLS



The Harris Uni-Core is a disposable multi-purpose sampling tool for soft substrates. Made from polypropylene, the cores have a hollow body with a built-in steel razor sharp cutting tip and a plunger for releasing the punched sample. Uni-Cores are individually wrapped and sterilized by ethylene oxide. Uni-Cores can be reused by rinsing the tip with ethanol after each use and autoclaving for 20 minutes at 250°F/15 psi.

Size:	Item#:
0.5mm	SIP0270
0.75mm	SIP0271
1.0mm	SIP0272
1.25mm	SIP0273
2.0mm	SIP0275
2.5mm	SIP0276
3.0mm	SIP0277

Size:	Item#:
3.5mm	SIP0278
4.0mm	SIP0279
5.0mm	SIP0280
6.0mm	SIP0281
7.0mm	SIP0282
8.0mm	SIP0283

HARRIS CUTTING MAT™

These versatile, inert, self-healing polyvinyl rubber cutting mats are designed for repeated, extended use and easy maintenance. Each mat offers two usable surfaces and the self-healing characteristics conceals cuts. The Harris Cutting Mat™ may be periodically cleaned with ethanol. The Harris Cutting Mat™ is available in two sizes and comes in a pack of five. Made for use with the Harris Uni-Core Tools.



Size:	Item#:
7.5x6cm	LWC0601P
20.5x15.5cm	LWC0602P



American MasterTech
scientific laboratory supplies

Copyright ©2011 American MasterTech Scientific Incorporated. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.