



American MasterTech
scientific laboratory supplies



HISTOPress

EMBEDDING MADE EASY



- Available in Pink, Purple, or White
- Xylene Resistant
- Lightweight and easy to handle with a cool shaft for holding
- Flattens tissue specimens evenly with uniform margins
- Fits small, medium, and large embedding molds

When a piece of embedded tissue is placed into a cassette mold that contains molten paraffin, it is very important to have the tissue as flat as possible so that when sectioning, the ribbons have even pieces of tissue. With the HistoPress tissue tampers, flattening the tissue within the cassette mold can be done easily. This great tool is designed with an innovative shaft (handle) that is long enough to hold comfortably while staying cool even with the base being hot. The HistoPress is a great addition to any histotech's embedding station.



Size:	HistoPress (Each)			HistoPress (Pack of 2)		
Description	White Item#	Pink Item#	Purple Item#	White Item#	Pink Item#	Purple Item#
Round (1/2")	TCZ1300	TCZ1300PK	TCZ1300PR	TCZ1300P	TCZ1300PKP	TCZ1300PRP
Small (1/2 x 1/2")	TCZ1302	TCZ1302PK	TCZ1302PR	TCZ1302P	TCZ1302PKP	TCZ1302PRP
Medium (3/4 x 9/16")	TCZ1304	TCZ1304PK	TCZ1304PR	TCZ1304P	TCZ1304PKP	TCZ1304PRP
Large (3/4 x 1-1/8 in")	TCZ1306	TCZ1306PK	TCZ1306PR	TCZ1306P	TCZ1306PKP	TCZ1306PRP

Copyright ©2011 American MasterTech Scientific Incorporated. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.